



The Accutome A-Scan Plus Connect® Software is DICOM Compatible.

Item Number: 24-4400

The DICOM (Digital Imaging and Communications in Medicine) is the standard for communication of medical imaging devices, by specifying a protocol for network communications, media communications, file formats and other information.

DICOM Capabilities

- Export and store to your DICOM server or PACS/ image management system.
- Verification of connection to other Application Entities (AEs).
- The only A-Scan to export axial length data along with reports to the DICOM/PACS, allowing manipulation of data and calculations within your DICOM system.

Features

- One license per facility address.
- Easy remote activation.
- Easy exporting.

As always, our product specialists are here to support your installation and training.

3222 Phoenixville Pike • Malvern PA 19355 USA • 800.979.2020
610.889.0200 • F:610.889.3233 • www.accutome.com • info@accutome.com

 **ACCUTOME**[®]
A HALMA COMPANY