




Welch Allyn® RetinaVue® Network

Secure transfer of encrypted retina images and management of retinal exam data with award-winning, HIPAA compliant RetinaVue Network software.

RetinaVue Network software offers a variety of PC- and web-based applications including comprehensive population health management and quality reporting tools to more effectively manage retinal exam data by clinic and by individual patient.



To streamline documentation, fully integrated, bi-directional interfaces with EMRs including Allscripts, athenahealth, Cerner, Epic, NextGen, and many others are offered. Place retina exam orders and automatically access diagnostic reports from your EMR.

RetinaVue Network Features

Web-based Portal Software



- Schedule a Patient (wireless only)
- View Scheduled Patients (wireless only)
- View New Diagnostic Reports
- Search Diagnostic Reports
- View Statistics
- Manage Devices (wireless only)
- Manage Clinics
- Edit Company Information

PC-Software Application



- Schedule a Patient
- View Scheduled Patients
- Transmit Retinal Images
- Manage General Settings such as exam state, clinic contact information, etc.
- Manage Camera and Network Settings
- Manage Contact Settings
- Edit Advanced Settings

A RetinaVue Network software subscription is required and priced per camera, per month.

RetinaVue Network System Requirements

Item	System Requirements
Configuration	Tower or desktop PC preferred; laptop or tablet
CPU	1 GHz (or equivalent) minimum, 2 GHz or greater recommended
Hard Disc	150 MB free HDD space, 16 GB free HDD space or greater recommended
RAM	2 GB minimum, 4 GB recommended
Monitor Resolution	1280 x 720 recommended
Ethernet Port	RJ-45
Ports	2 USB, 2.0 port or greater
Operating System	Windows® 7 with SP1 32-bit and 64-bit, Windows 8.1 64-bit, Windows 10 with latest SP RetinaVue Network Prerequisite software requirements: <ul style="list-style-type: none"> • Microsoft Visual C++ 2013 Runtime Libraries (x86) • Microsoft® .NET Framework 4.5 • Perform a Windows update before installation
Web browser	Internet Explorer® (version 11 or greater), Chrome™ for Windows (latest version)
High-Speed Internet Connection	Broad-band Internet connection (minimum download speed 1.5 Mbps)
Firewall	Ability to connect to www.retinavue.net on ports 80 and 443 with RetinaVue Network.exe <ul style="list-style-type: none"> • Allows the submittal of captured images from the client to the Welch Allyn RetinaVue Network Web Server
Certificates*	Go Daddy® Root Certificate Authority installed on the computer <ul style="list-style-type: none"> • Allows secure/protected communication between the client and the Welch Allyn RetinaVue Network Web Server • This is usually installed by default on the supported Windows Operating system for this product. Contact your IT department for questions and issues.

*IT Network Security: Computers used to connect to RetinaVue Network Customer Portal or used to download the RetinaVue Network software application should be set up and maintained following IEC 80001 or similar IT Network security practices.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361