

No Slip Trial Frame (PAT.) by Hasegawa BicoH

MD674(PD54~74)

NEW

Max: P D 74



No slip temple tip. Max PD size 74mm.

- Please move the dial. You are able to change PD easily.
- 11 sizes (PD 54 ~ 74)
- No Slip Temple Tip (PAT.)
- Light! Less 50g
- You can select the other type as follows..



straight temple



half rim(for progressive lens)

MD547(PD54~70)

Renewal



We end sales as soon as this stock ends.

- Please move the dial. You are able to change PD easily.
- 9 sizes (PD 54 ~ 74)
- No Slip Temple Tip (PAT.)
- Light! Less 50g
- You can select the other type as follows..



straight temple



half rim(for progressive lens)

MD486(PD48~60)

Renewal



- Please move the dial. You are able to change PD easily.
- 7 sizes (PD 48 ~ 60)
- No Slip Temple Tip (PAT.)
- Light! Less 50g
- You can select the other type as follows..



straight temple



half rim(for progressive lens)

No slip temple tip

NEW



- It is possible to change for your 486 and 547.
- When you exchange it for this temple tip, slipping is reduced.