

Item Number: 24-6100

The B-Scan Plus' proprietary software combined with unique probe electronics make it the only portable High Definition B-Scan. While developing the B-Scan Plus software, Accutome visited many of the top facilities around the world. The results are a reliable and easy-to-use unit, which can quickly scan patients and transport information.

Features:

- Image is Everything 0.015 mm resolution is the highest in the medical industry.
- State-of-the-Art Probe Design Sharper, more focused images due to the elimination of signal loss.
- Unsurpassed Data Analysis Contains 6 useful measuring calipers.

- "Smooth Zoom" Technology The only unit with 2x full image zoom without distortion of real-time or captured scan.
- Portable Plug probe into any Windows based laptop or PC and go.
- Never Miss an Image Unlimited 34 second film loops.
- Share Information Easily Adaptable document transfer via EMR, email or printer.
- Compact File Storage Unsurpassed by the competition.
- Protect Your Investment Easily upgradable software.
- Master in Minutes User-friendly interface makes user manuals and quick start guides unnecessary.

